



MEDICALMART®

Verify® Self-Contained Biological Indicators (SCBIs) are for use in monitoring steam sterilization and ethylene oxide (EO) sterilization processes. Chemical indicators on each vial provide immediate proof of processing. Following incubation, a VIVID COLOUR CHANGE from deep blue/purple to bright yellow and/or turbidity gives unmistakable evidence of microbial growth. 24-hour final results for steam sterilization & 48-hour final results for EO sterilization

STERIS® VERIFY® DUAL SPECIES SCBI BIOLOGICAL INDICATORS

- 779-S32716 Well Incubator Each
- 779-S3075Activator Each
- 779-S3061Biological Indicators Box/50

24 HR Readout



FREE INCUBATOR WITH THE
PURCHASE OF 3 BOXES OF
VERIFY® SCBIs

STERIS® VERIFY® STEAM INTEGRATING INDICATOR STRIPS, CLASS 5

The Verify Integrating Indicators are used to monitor steam sterilization within packs. Each indicator consists of a plastic carrier with a moving front indicator to illustrate a pass or fail condition.

- Internal pack monitoring
- 250°F/121°C 30-minute gravity
- 270°F/132°C 3-minute gravity IUSS
- 270°F/132°C 4-minute prevacuum

- 779-PCC067 Indicators for steam, 2" Box/500
- 779-PCC064 Indicators for steam, 4" Box/250
- 779-PCC064B Indicators for steam, 4" Case/2000



779-PCC067



779-PCC064

STERILIZATION POUCHES, SELF SEALING

New regulations require sterilization pouches to be Class 4 (Highest Level) which these pouches meet. Requires a Class 5 indicator strip (below) to meet new regulatory guidelines for processing.

- 211-SCXX2 2.25" x 4" Box/200
- 211-SCX2 2.25" x 9" Box/200
- 211-SCXS2 3.5" x 5.25" Box/200
- 211-SCS2 3.5" x 9" Box/200
- 211-SCM2 5.25" x 10" Box/200
- 211-SCL2 7.5" x 13" Box/200



Ask us for a Free Consultation on your Sterilization Procedures today!

www.medicmart.com

Ontario & Quebec Phone (905) 624-6200 1-800-268-2848
Atlantic Canada Phone (902) 468-4087 1-800-268-2848

Newfoundland Phone (709) 726-5877 1-888-842-9462
Western Canada Phone (780) 444-5800 1-866-444-5848